



COBI REHAB®
The Positioning Company

XXL-Rehab Stok

XXL-Rehab Walking Cane

Anvendelse
Use

Max 270 kg

Brugsanvisning
User manual

DK

Kan med fordel anvendes til:

- Bariatriske brugere.

UK

Ideal use for:

- Bariatric users.



ENDynamics Ltd.
1 Gladstone Court
Farnworth, Bolton, BL4 7EU

Tel: 0333 444 1034
Email: email@endynamics.co.uk

Maintenance and cleaning Specifications

The XXL-Rehab Walking Cane can be cleaned with water and ordinary detergent.

Clip and the continuous rod for height adjustment should be regularly inspected and if necessary replaced.

⚠ Failure to inspect the clip may cause malfunction.

Maximum weight capacity is	270 kg
Height adjustable	71-93 cm
Weight	0.9 kg
Material	Anodized steel
Grip	Vinyl
ISO - classification	ISO 12 03 03 - 05

Warranty

The XXL-Rehab Walking Cane is covered by one year warranty from the invoice date for materials and workmanship.

Products which have been used for purposes other than the intended use; or which have been subject to negligence, abuse, improper storage or handling, improper operation, unauthorized modifications or damages beyond normal wear and tear as determined by Cobi Rehab, are not covered by this warranty.

If weight capacity is exceeded, the warranty will be void.

Any unauthorized repairs to product / part, as well as tampering with any components, will void the warranty.



XXL-Rehab Walking Cane is manufactured and tested in accordance with international standards EN ISO 11334-1:2007.



MANUFACTURER AND DEALER

Cobi Rehab
Fuglebækvej 1D
DK-2770 Kastrup
Denmark
Tel.: +45 70 25 25 22
www.cobirehab.com

04/2018/JP